

Evolution of the AffloVest



1st Generation
(2013)



(2014)



(2016)



(2017)



5th Generation
(2020)

Original
Controller &
Travel Bag



Controller &
Travel Bag
Since 2016



Continuous Innovation Has Driven Improved Product Performance & Market Leading Technology
Providing more innovation in the past 7 years than the 30-year history of HFCWO

- ✓ **Strength:** vibratory force doubled since the original AffloVest – ability to adjust power settings to meet any patient’s needs
- ✓ **Controller:** upgraded from original analog controller to current digital controller with ability to track compliance
- ✓ **Ergonomics:** custom fit and 28% weight reduction – AffloVest is the lightest, best-fitting HFCWO vest in the market
- ✓ **Battery Life:** consistent improvements in battery life – current AffloVest battery can last ~2 hours on the high setting (or ~2-4 days of use⁽¹⁾)

1. Dependent upon frequency of patient use – range of 2 - 4 days assumes 30 - 60 minutes of use per day

“My patient can tolerate the highest frequencies of all three settings with the AffloVest feeling more intense on her lungs than her current compression system...”
– Pediatric Pulmonologist at a Children’s Hospital

Feel the Difference
For more information, please visit [AffloVest.com](https://www.afflovest.com)